

magnetiQ

Total RNA Extraction Kit

Quick Start Guide

1 In 1.5ml microfuge tube, add 380µl of **Lysis/Binding Buffer** and then add 20µl of binding magnetic nanoparticles.




2 Add 100µl of cell suspension in PBS containing a maximum 10^7 cells to the microfuge tube.

Mix cells and buffer well by pipetting up-down 10x, then **incubate 2 mins** to allow for lysis and RNA binding.




3 Place tube on magnetic rack for 2 mins to capture the RNA-bead complex, then discard supernatant.



4 Remove tube from magnetic rack and **resuspend RNA-bead complex in 600µl of Wash Buffer**. Mix well by pipetting up-down 10–15x.

Return to magnetic rack for 1 min, then discard supernatant.



5 Repeat wash with 600µl of Wash Buffer.

Return to magnetic rack for 1–2 mins, discard supernatant and leave to dry for 2 min.




6 Remove tube from magnetic rack and **resuspend RNA-bead complex in 50µl of Elution Buffer**, then mix well by pipetting up-down 15–20x to elute RNA from beads.

Let stand for an additional 1 min.



7 Place tube on magnetic rack to separate beads (~2 mins).



8 Transfer clean RNA solution to clean tube.

